GX-L/GF-L

Quick Start Guide GX-L/GF-L Series High Capacity Precision Balance	English (EN)	3
Guide de Démarrage Rapide Série GX-L/GF-L Balance de Précision de Grande Capacité	Français (FR)	14
Kurzanleitung Baureihe GX-L/GF-L Hochleistungs-Präzisionswaage	Deutsch (DE)	26
Guía de Inicio Rápido Serie GX-L/GF-L Báscula de Precisión de Alta Capacidad	Español (ES)	38
Guida Introduttiva Rapida Serie GX-L/GF-L Bilancia di Precisione ad Alta Capacità	Italiano (IT)	50
快速入門指南 GX-L/GF-L 系列 高秤量精密天平	中文 (CT)	62

A&D Company, Ltd.

Note for Precautions

 **CAUTION** A potentially hazardous situation which, if not avoided, may result in personal injury or property damage.

CAUTION Describes the points to be careful for appropriate use.

NOTE Describes 'highly possible to be handled inappropriately' or 'general advice in using the device'

About This Manual

- (1) No parts of this manual may be reproduced without permission. This manual may not be copied, modified, or translated without the written permission of A&D Company, Limited.
- (2) The contents of this manual are subject to change without notice.
- (3) Please contact A&D if you notice any uncertainty, errors, omissions, etc. in this manual.
- (4) A&D bears no liability for direct, indirect, special, or consequential damages due to the operation of this product, even if advised of the possibility of such damage. Furthermore, A&D assumes no liability for claims of rights from third parties. Concurrently, A&D assumes no liability whatsoever for data losses regardless of (3) above.

1. Introduction

Thank you for purchasing a GX-L/GF-L series precision balance from A&D. This manual describes the basic functions and operations of the balance. For more details, refer to the manual described in "1-1. List of Functions Explained in a Separate Manual".

1-1. List of Functions Explained in a Separate Manual

The functions and operating procedures of the balance are explained in detail in documents that supplement this manual. Download the following from our website, <https://www.aandd.jp>.

[GX-L/GF-L Series Instruction Manual]

- Instruction manual for fully understanding and utilizing GX-L/GF-L series precision electronic balances.


Scan this QR code for the manual.

2. Installation and Precautions

The balance is a precision instrument. Unpack it carefully. The packing contents differ by balance model. Refer to the images on the following page as you check to make sure all items are present. Keep the packing materials to use when transporting the balance in the future.

(1) Standard type

GX-12001L/22001L/32001L/32001LD/42001L/62001L/62000L/102000L
 GF-12001L/22001L/32001L/62000L

NOTE


- When moving the arm, loosen the knob. When fixing the arm, fasten the knob.

⚠ CAUTION

- Use the dedicated AC adapter for the balance.
- Do not use the included AC adapter for other devices.
- Using the incorrect AC adapter may cause performance issues in the balance or other devices.
- The dustproof and waterproof structure may cause difficulty when inserting the USB cable, so make sure it is fully plugged in. Put on a cap when not in use.
- The dustproof and waterproof structure may cause difficulty when inserting the AC adapter, so make sure it is fully plugged in.
- When using the balance under conditions requiring dustproof and waterproof performance, attach the terminal cover or use a waterproof RS-232C cable (AX-KO2737-500).

(2) Removable Display Unit Type

GX-32001LS/32001LDS/62001LS/62000LS/102000LS


⚠ CAUTION

- ❑ Use the dedicated AC adapter for the balance.
- ❑ Do not use the included AC adapter for other devices.
- ❑ Using the incorrect AC adapter may cause performance issues in the balance or other devices.
- ❑ If you have multiple balances, make sure the serial numbers of the balances and display you intend to use are identical before connecting them via the cable. Failure to do so may result in performance issues.
- ❑ The dustproof and waterproof structure may cause difficulty when inserting the USB cable, so make sure it is fully plugged in. Put on a cap when not in use.
- ❑ The dustproof and waterproof structure may cause difficulty when inserting the AC adapter, so make sure it is fully plugged in.
- ❑ When using the balance under conditions requiring dustproof and waterproof performance, attach a terminal cover or use a waterproof RS-232C cable (AX-KO2737-500).

2-1. Installation Considerations, Preparation and Precautions

Prepare the following installation conditions in order to ensure optimal balance performance.

- ❑ The best operating temperature is about $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$ at about 45% to 60% relative humidity.
- ❑ Install the balance in a location that is free of dust.
- ❑ Install the balance on a firm and horizontal floor (not a soft floor).
- ❑ The weighing table should be solid (an anti-vibration table or stone table is ideal).
- ❑ Make sure that the location where the balance is installed is not at an angle.
- ❑ Avoid vibration. Corners of rooms on the first floor are best, as they are less prone to vibration.
- ❑ Install the balance where it is not affected by breezes from air conditioners, etc.
- ❑ Install the balance where it is not exposed to direct sunlight.
- ❑ Install the balance away from equipment which produces magnetic fields.
- ❑ Level the balance by adjusting the leveling feet so that the bubble is centered in the black circle of the bubble spirit level.
- ❑ Make sure to warm up the balance before use for at least an hour with the AC adapter connected.
- ❑ To ensure weighing accuracy, adjust the sensitivity of the balance with the AC adapter connected before using it for the first time or after moving it to another location. For how to adjust sensitivity, refer to "7. Sensitivity Adjustment/Calibration Test" in the separate instruction manual.
- ❑ The dustproof and waterproof level of the balance is equivalent to IP65, and its second digit, "5", corresponds to "having no harmful influence from direct jets of water." Washing the balance with high-pressure water or submerging it in water may cause water to enter the balance, leading to malfunction.
- ❑ When installing and using the balance under conditions requiring dustproof and waterproof performance, make sure that the AC adapter plug is fully inserted into the AC adapter jack and that the terminal cover is attached to the RS-232C interface or use a waterproof RS-232C cable (AX-KO2737-500).
- ❑ If the RS-232C terminal cover is removed or you are not using a waterproof RS-232C cable (AX-KO2737-500), the balance will not be protected from dust or water.


⚠ CAUTION


Do not install the balance where corrosive or flammable gas is present.

3. Display Symbols and Key Operations

Display symbols


Blinking display contents


No.	Description
1	Processing indicator
2	Stabilization indicator
3	USB connecting indicator
4	Lit: Power source on standby Blinking: Automatic sensitivity adjustment
5	Displays weighing value, stored data, and item name
6	Units
7	Net indicator
8	Gross indicator
9	Number of static data (statistical calculation mode) Data memory number Displays the percentage of weight data relative to the weighing capacity (capacity indicator) Displays the setting value of the function settings
10	Response indicators (light for approximately 30 seconds after weighing starts)
11	Comparator indicators
12	Display holding indicator
13	ISD shock indicator
14	Gross zero indicator
15	Lit: Interval memory on standby Blinking: Interval memory active
16	Preset tare indicator

Key operations

Key operations affect how the balance functions. Normal key operations during measurement are "press and release the key immediately" or "press the key". Do not "press and hold the key (for two seconds) unless specifically required.

Press the key
(Press and release
the key immediately)

Press and hold the key
for 2 seconds

Key	When pressed and released	When pressed and held (for 2 seconds)
	Turns the display on and off. Only the standby indicator is displayed when the display is turned off. Weighing mode is enabled when the display is turned on. When the password function is enabled, the password input display is displayed. Refer to "16-2. How to Input the Password at the Start of Weighing" in the separate instruction manual. This ON:OFF key is always enabled. Pressing the ON:OFF key during operation will interrupt operation and turn the display off.	
	Adjusts the sensitivity of the balance using the internal mass (GX-L series only).	Displays other sensitivity adjustment menu items.
	Switches the weighing unit stored in the function table.	Enters the self check function mode. Refer to "6-2. Self Check Function/Automatic Setting of Minimum Weighing Value by ECL".
	<ul style="list-style-type: none"> <input type="checkbox"/> In weighing mode, turns the minimum weighing value on and off. <input type="checkbox"/> In counting or percent mode, enters sample storing mode. 	<ul style="list-style-type: none"> <input type="checkbox"/> Enters function table mode. Refer to "9. Function Table". <input type="checkbox"/> Runs the repeatability check function when pressed and held for another 2 seconds after the function table menu is displayed (GX-L series only). Refer to "17. Repeatability Check Function" in the separate instruction manual.
	Depending on function table settings, either stores the weighing data or outputs the data to a printer when the value is stable (factory default = output).	<ul style="list-style-type: none"> <input type="checkbox"/> Enters the mode for changing the unit mass registration number in counting mode. <input type="checkbox"/> Depending on changes to the function table: <ul style="list-style-type: none"> • Outputs "Title block" and "End block" for GLP/GMP reports. • Displays the data memory menu. • Enters the mode for the density number in flow measurement.
	Sets the display to zero. Refer to "About zero-point operation", "About tare", "About measurement range" in "4-2. Basic Operation (g, kg, ct, mom)" in the separate manual.	
		

4. Sensitivity Adjustment

4-1. Automatic Sensitivity Adjustment (GX-L Series Only)

This function automatically adjusts the sensitivity of the balance according to ambient temperature change, set time or interval time using the internal mass. For details, refer to "7-1. Automatic Sensitivity Adjustment" in the separate manual.

4-2. Sensitivity Adjustment Using the Internal Mass (GX-L Series Only)


Sensitivity adjustment using the internal mass can be performed with one key press.

1. Connect the AC adapter and warm up the balance for at least 30 minutes with nothing on the weighing pan.
2. Press the **[CAL]** key. The balance displays **[CAL in]**.
3. The sensitivity of the balance is adjusted using the internal mass. Do not allow vibration or drafts to affect the balance.
4. After sensitivity adjustment, if GLP output is set, a "sensitivity adjustment report" is output or stored to data memory.
5. The balance returns to weighing mode automatically.

About the internal mass

The value of the internal mass may change due to factors such as the operating environment and aging. Correct the internal mass value as necessary. Refer to "7-6 Correcting the Internal Mass Value (GX-L Series Only)" in the separate manual.

Since the internal mass is about 1 kg, the possibility of error may increase as the weighing value increases. To maintain the weighing accuracy, perform the sensitivity adjustment using an external weight periodically, as described in "4-3. Sensitivity Adjustment Using an External Weight".


This indicates that the balance is importing calibration data. Do not allow vibration or drafts to affect the balance while this indicator is displayed. After calibration, the balance returns to indicate the previous display.













4-3. Sensitivity Adjustment Using an External Weight














This function adjusts the sensitivity of the balance using an external weight.

1. Connect the AC adapter and warm up the balance for at least 30 minutes with nothing on the weighing pan.
 2. Press and hold the **[CAL]** key for 2 seconds.
 3. Continue holding the key until **[CAL out]** is displayed. Once **[CAL out]** is displayed, release the key.
 4. Make sure that nothing is on the weighing pan and press the **[PRINT]** key.
 5. The balance measures the zero point. Do not allow vibration and the like to affect the balance.
 6. The balance displays the weight value (g).
 7. Place an external weight on the weighing pan and press the **[PRINT]** key. The balance measures the external weight.
 8. The balance measures the sensitivity adjustment value. Do not allow vibration and the like to affect the balance.
 9. Remove the external weight from the weighing pan.
 10. After sensitivity adjustment, if GLP output is set, a "sensitivity adjustment report" is output or stored to data memory.
 11. The balance returns to weighing mode automatically.
 12. Place the external weight on the weighing pan again to check if it is within the set value ± 2 digits^{*1}. If it is not within that range, start over from the first step of this procedure in the appropriate ambient conditions.
- *1 "Digit" is the minimum display unit.
For example, if minimum display is 0.1 g, 1 digit is equal to 0.1 g.

5. Error Codes

Display	Error code	Description
		Internal error If this error continues to be displayed, please request repair from your local dealer.
	EC, E11	Stability error The balance is unable to stabilize due to the environment. Check around the pan. Prevent vibration, drafts, temperature changes, static electricity and magnetic fields from affecting the balance. Refer to "2-3. During Use" in the separate manual. To return to weighing mode, press the  key.
		Outside of the setting range The data to be stored is outside of the setting range. Set again within the setting range.
		Balance internal memory element malfunction If this error continues to be displayed, please request repair from your local dealer.
	EC, E16	Internal mass error (GX-L series only) Applying the internal mass does not yield a change in the mass value as specified. Confirm that there is nothing on the pan and perform the operation from the beginning again.
	EC, E17	Internal mass error (GX-L series only) The internal mass application mechanism does not function properly. Perform the operation from the beginning again. If the problem persists, please request repair from your local dealer.
		Abnormality in the internal memory data of the balance If this error continues to be displayed, please ask for repair to your local dealer.
		Balance internal memory abnormality If this error continues to be displayed, please request repair from your local dealer.
	EC, E20	Sensitivity adjustment weight error (Positive value) The sensitivity adjustment weight is too heavy. Check around the pan. Check the mass value of the weight. To return to the weighing mode, press the  key.
	EC, E21	Sensitivity adjustment weight error (Negative value) The sensitivity adjustment weight is too light. Check around the pan. Check the mass value of the weight. To return to the weighing mode, press the  key.
		Overload error A sample beyond the balance weighing capacity has been placed on the pan. Remove the sample from the pan. If the problem persists, please request repair from your local dealer.
		Weighing pan error The mass value is too light. Check that the weighing pan is installed correctly. Set the weighing pan correctly. Adjust the sensitivity of the balance.
		Overload error A sample beyond the balance weighing capacity has been placed on the pan. Remove the sample from the pan. If the problem persists, please request repair from your local dealer.
		Weighing pan error The mass value is too light. Check that the weighing pan is installed correctly. Set the weighing pan correctly. Adjust the sensitivity of the balance.

Display	Error code	Description
		Sample mass error The sample is too light to be stored as a sample mass for counting mode or percent mode. Use a larger sample.
  		Unit mass error The sample unit mass for counting mode is too light. Storing and using it for counting will cause a counting error. Add samples to reach the number indicated and press the PRINT key. Pressing the PRINT key without adding samples will put the balance in counting mode, but to ensure accurate counting, be sure to add samples.
		Clock battery error The clock backup battery has been depleted. Press any key and set the time and date. Even if the clock backup battery is depleted, the clock and calendar function works normally as long as the balance is powered with the AC adapter. If this error appears frequently, contact your local A&D dealer.
		Power supply voltage fault The voltage supplied from the AC adapter is abnormal. Make sure you are using the AC adapter supplied with the balance.
 		ECL repeatability error With the self-check function, the standard deviation (SD) of repeatability with electronically controlled load (ECL) exceeded 50 digits*1. Check the installation environment of the balance. <input type="checkbox"/>  ECL repeatability error <input type="checkbox"/>  ECL minimum weighing value (reference value) error Refer to "6-2. Self Check Function/Automatic Setting of Minimum Weighing Value by ELC" in the separate manual.
 ↑ Alternate (Blink) 		Full memory The number of stored weighing values has reached the upper limit. In order to store new weighing values, it is necessary to delete data. Refer to "11. Data Memory" in the separate manual.
 ↑ Alternate (Blink) 		Full memory The number of items stored in your sensitivity adjustment/calibration history has reached 50. When new items are added, older items are deleted. Refer to "11. Data Memory" in the separate manual.
		Level check Check the bubble level on the main unit. Rotate the leveling feet to center the bubble in the circle. Either adjusting the bubble level or pressing any key other than the ON:OFF key will clear the warning.

*1 "Digit" is the minimum display unit.

For example, if minimum display is 0.1 g, 1 digit is equal to 0.1 g (50 digits are 5.0 g).

6. Specifications

GX-L series (with internal mass)	GX-12001L	GX-22001L	GX-32001L GX-32001LS	GX-32001LD GX-32001LDS	GX-42001L
GF-L series (general)	GF-12001L	GF-22001L	GF-32001L		
Weighing capacity	12 kg	22 kg	32 kg	32 kg	42 kg
Maximum display	12.0084 kg	22.0084 kg	32.0084 kg	32.008 kg 6.2009 kg *1	42.0080 kg
Readability	0.1 g			1 g / 0.1g	0.5 g
Repeatability (Standard deviation)	0.1 g			0.5 g / 0.1 g	0.5 g
Linearity	±0.2 g			±1 g / ±0.2 g	±1 g
Accuracy right after sensitivity adjustment using the internal mass*2	±1.0 g		±1.5 g		
Dust and water protection	Complies with IP65				

GX-L series (with internal mass)	GX-62001L GX-62001LS	GX-62000L GX-62000LS	GX-102000L GX-102000LS
GF-L series (general)		GF-62000L	
Weighing capacity	62 kg	62 kg	102 kg
Maximum display	62.0084 kg	62.084 kg	102.084 kg
Readability	0.1 g	1 g	1 g
Repeatability (Standard deviation)	0.2 g	0.7 g	1 g
Linearity	±0.5 g	±1 g	±2 g
Accuracy right after sensitivity adjustment using the internal mass*2	±3.0 g	±5 g	±10 g
Dust and water protection	Complies with IP65		

*1 Weighing is enabled in a precision range, even with high tare values (smart range function).

*2 Values when there are no adverse effects due to abrupt changes in temperature, drafts, magnetic fields or static electricity (GX-L series only).

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A&D AUSTRALASIA PTY LTD
A&D KOREA Limited
A&D Instruments India Private Limited