

DIGITAL BLOOD PRESSURE MONITOR

UA-767 Plus

Clinically
Validated

IHB



The Monitor Tells You More

- Irregular Heart Beat (IHB) Indicator
- Blood Pressure Classification Indicator
- Last Reading Recall
- Extra Large 3-line Display for an Easier Read
- One-touch Measurement
- SlimFit Comfortable Cuff

CE 0123

A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

DIGITAL BLOOD PRESSURE MONITOR UA-767 Plus

Irregular Heart Beat (IHB) Indicator

Lights up after IHB is detected during measurement. Informs users of an irregular heartbeat, which may cause cardiovascular disease.

(Note: This is for warning purposes only. A&D recommends contacting a physician if the user sees this indicator light frequently)

Blood Pressure Classification Indicator

World Health Organization (WHO) classifies blood pressure ranges into 6 stages. UA-767 Plus indicates your position (Optimal to Severe hypertension) automatically after each measurement is complete.

Last Reading Recall

UA-767 Plus indicates your previous reading on the LCD display when the START button is pressed.

Extra Large 3-line Display for an Easier Read

Easy to read 3-line display shows all the measurement (systolic/diastolic pressure and pulse rate) simultaneously. Larger numbers are easier on the eyes.

One-touch Measurement


Pressing the START button completes all the steps for a blood pressure measurement, for simple and accurate operation.

SlimFit Comfortable Cuff

A new innovative edge free internal bladder "EFB" promises even pressure on the arm during the measurement. This provides the user with a higher level of comfort while the arm is under pressure.

WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart below.



A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

*UA767PlusCA-ADCC-14-PR6-143**
REF*CARD-UA767Plus