

DIGITAL BLOOD PRESSURE MONITOR

UB-512

Clinically
Validated

IHB



Family Wrist Monitor

- Ultra Compact and Light Weight Design
- Fast Measurement along with Inflation
- Irregular Heart Beat (IHB) Indicator
- Blood Pressure Classification Indicator
- Twin 30 Memory + Average Reading
- Illuminated Start Button

CE 0123

A&D
A&D Medical

Technology You Can Trust

<http://www.aandd.jp>

Ultra Compact and Light Weight Design

Our state-of-the-art technology brings real portability to the wrist monitor. It is 40% smaller and 30% lighter than our former model. Such miniature monitor promises portable measurement whenever user needs to monitor blood pressure.

Fast Measurement Along with Inflation

The UB-511/512 series completes blood pressure measurement while inflating the wrist cuff. This gives user a quick measurement at anywhere he needs to monitor his blood pressure.

Irregular Heart Beat (IHB) Indicator

This lights up if IHB is detected during the measurement. Helps users to find an irregular heartbeat which may be associated with cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user sees this indicator frequently)

Blood Pressure Classification Indicator

World Health Organization (WHO) classifies blood pressure ranges into 6 stages. This monitor indicates your position (Optimal to Severe Hypertension) automatically after each measurement is complete.

Twin 30 Memory + Average Reading

The UB-512 has two sets of large memory for 30 readings for monitoring the trend of measurements. An average of all the stored data is shown at each start up, giving a condensed reading history at a glance.

Illuminated Start Button

A&D has included an innovative idea for ease-of-use for two people sharing the monitor. The UB-512 has illuminated buttons and a button (1 or 2) lights up while measurement is in progress.

WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), as shown in the chart on the right.

