

# Bedside Scale

Model: AD-6121A

The AD-6121A is a lightweight, portable scale ideal for bedside use.

## General Specifications

- Can display measurements in kg
- Lightweight 3 kg design with handle for portability and convenient storage
- Use at the bedside and store underneath
- Adjustable indicator
- Auto power off function
- Supports AA battery operation or AC adapter connection



RS-232C  
Interface



LCD  
Display



Auto  
Power Off



Normal  
Battery



Lightweight & compact size



Level adjustment



## Technical Specifications

Capacity	150 kg
Graduation	0.1 kg
Measurement Range	2 kg to 150 kg
Power Source	4 Batteries LR6 (AA) / AC Adapter
Battery Life	Approx. 1,000 measurements (Battery use)
Interface	RS 232C 1ch (3Pin)
Dimensions (W) x (D) x (H)	315 x 314 x 57 mm
Weight	3 kg (Unit without batteries)
Functions and Properties	TARE, Auto power-off

## Optional Accessories

Product Code	Option Name
AX-KO2238	Communication cable
AD-6121-01	Carry case



Communication  
cable



Carry case

## What is In The Box

- 1 Instruction manual
- 1 Warranty

- Please purchase an AC adapter that complies with the standards of your region.
- Ensure that the product complies with the regulatory requirements of the country where it will be used before use.
- Specifications may differ in each country. Please contact an A&D subsidiary for more information.

# Baby Scale “Baby One”

Model: AD-6020-12K

The AD-6020-12K is a baby scale with a removable basket for easy cleaning.

## General Specifications

- Can display measurements in kg
- Ensure hygiene with easy-to-clean removable basket
- Large 25 mm digital display
- Supports D battery operation or AC adapter connection
- Record and recall last measurement values with a single button
- Tare function



**Memory**  
Memory Function

**888**  
LCD Display


**OFF**  
Auto Power Off

**B+**  
Normal Battery

## Height measurement

The basket has a scale and a cursor for easy height measurement. The cursor for height measurement can be stored in the lower part of the main unit.

Cursor for height measurement



## Technical Specifications

Capacity	12 kg
Graduation	10 g < 10 kg > 20 g
Measurement Range	0.1 kg to 12 kg
Height Measurement Range	45 cm to 57 cm
Power Source	4 Batteries LR20 (D) / AC Adapter
Dimensions (W) x (D) x (H)	628 x 343 x 245 mm
Weight	5 kg (Unit without batteries)
Functions and Properties	TARE, Height measurement, Safety lock



## What is In The Box

- 1 Height measurement cursor
- 1 Instruction manual
- 1 Warranty

- Please purchase an AC adapter that complies with the standards of your region.
- Ensure that the product complies with the regulatory requirements of the country where it will be used before use.
- Specifications may differ in each country. Please contact an A&D subsidiary for more information.