

Barrier-free Scale (Scale with Handrail)

Model: AD-6107R

The AD-6107R is a scale with a handrail to support patients during weighing. It is brightly colored to ensure correct positioning and is equipped with large, lightweight casters for portability.

General Specifications

- Can display measurements in kg
- RS-232C expansion terminal
- Battery operation for over 200 hours
- BMI display function
- Tare function

Easy to carry

- Large casters for easy mobility
- Long-life battery operation for quick measurement



Safety Design

- Easy-to-ride weighing platform with low 48 mm height
- Easy-to-read color identification of the weighing position
- Rounded corners for safety



Large handle for easy grabbing

The large handle allows you to hold on to the scale while measuring, ensuring accurate and safe weight measurement.



Technical Specifications

| | |
|--------------------------------|--|
| Approval Class | III (3) |
| Capacity | 200 kg |
| Graduation | 0.1 kg < 150 kg > 0.2 kg |
| Measurement Range | 2 kg to 200 kg |
| Power Source | AC Adapter(AX-TB294) / 4 Batteries LR20(D) |
| Interface | RC232C 2ch |
| Operating Conditions | 5 °C to 35 °C / Not more than 85 %RH (no condensing) |
| Transport / Storage Conditions | -10 °C to +60 °C |
| Dimensions (W) x (D) x (H) | 613 x 653 x 870 mm Platform: 493 x 426 x 48 mm |
| Weight | 17 kg (Unit without batteries) |
| Functions and Properties | TARE, BMI measurement |

Optional Accessories

| Product Code | Option Name |
|--------------|-----------------------------|
| AX-PP147-S | Printer paper(5 pieces/box) |
| AD-PR580 | AD printer unit |
| AD-PR580-007 | Printer installation kit |
| AS-BLS-120W* | AC Adapter for AD-PR580 |
| AS-ACS-120A* | AC code for AX-BLS-120W |

*This is an accessory for the AD-PR580.



Printer paper



AD printer unit

AD Printer Unit (AD-PR580)

Printer installation kit is required for mounting on the stand



Technical Data

| | |
|------------------|--------------------------------|
| Designated paper | AX-PP147-S W 58 mm x L 29 m |
| Power supply | AC adapter(AS-BLS-120W) |

What is In The Box

- 1 Sticker for AC adapter identification
- 1 Instruction manual
- 1 Warranty

- Please purchase an AC adapter that complies with the standards of your region.
- Ensure that the product complies with the regulatory requirements of the country where it will be used before use.
- Specifications may differ in each country. Please contact an A&D subsidiary for more information.



A&D Company Limited.(JAPAN)
URL:aandd.jp

A&D ENGINEERING, INC.(USA)
URL:andonline.com

A&D Instruments Canada Inc.(CANADA)
URL:autocontrol.com

A&D INSTRUMENTS LIMITED(United Kingdom)
URL:andprecision.com

A&D AUSTRALASIA PTY LTD(Australia)
URL:andaustralasia.com.au

A&D RUS CO., LTD.(Russia)
URL:and-rus.ru

A&D Technology Trading(Shanghai)Co., Ltd(China)
URL:aanddtech.cn

A&D INSTRUMENTS INDIA PRIVATE LIMITED(India)
URL:aanddindia.in